Product: **STENSIL ECO 70**

Knetsilikon / putty silicone

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Knetsilikon / putty silicone

Instruction for use

Please read carefully and retain for future reference

Working
room temperature
at 23° C

Delivery form:
2x1,43 Kg tins
2x4,77 Kg buckets

2x25 Kg buckets

Indications for use

STENSIL ECO 70 is an addition-curing 1:1 kneadable precision impression material with a final hardness of Shore A 70, its mixing ratio is very simple and save. STENSIL offers an extraordinary precise reproduction as a result of its detail graphic accuracy and a superior dimension stability due to a very low shrinkage, it is equipped with outstanding physical properties, high heat resistant, not sticky and physiologically uncritical.

STENSIL ECO 70 is perfect suitable for check bites, silicone keys, bite registrations, models for denture repairs, blocking-out, teeth-insulating in the flask and all kinds of moulding.

Working instructions

Take an equal part from both components (A white + B red) and knead it intensively and in a timely manner approx. 30 sec. till a strip-free and homogeneous mixture is produced, recognizably by an evenly light blue colour. Thereupon the mass can be appliquéd and after another approx. 3,5 min. it can be replaced.

Important working hints

- For with-drawal and mixing gloves made from vinyl not latex should be used since latex extends the setting time or prevents it possibly. If mixing bare-handed the setting time will be poss. cut down.

Technical datas

Dosing: 1:1
Density: 1,56 g / cm³
Colour code: blue
Hardness: Shore A 70
Mixing time: ca. 30 sec
Working time: ca. 1,20 minutes
Setting time: ca. 4,50 minutes
Linear dimension change: < 0,2%